

Development of reference materials to standardize microvesicle detection

Edwin van der Pol^{1,2}

on behalf of Rienk Nieuwland¹



June 29th, 2012

¹Laboratory Experimental Clinical Chemistry; ²Biomedical Engineering and Physics, Academic Medical Center, Amsterdam, The Netherlands



Disclosures for Edwin van der Pol


In compliance with COI policy, ISTH requires the following disclosures to the session audience:

Research Support/P.I.	No relevant conflicts of interest to declare
Employee	No relevant conflicts of interest to declare
Consultant	No relevant conflicts of interest to declare
Major Stockholder	No relevant conflicts of interest to declare
Speakers Bureau	No relevant conflicts of interest to declare
Honoraria	No relevant conflicts of interest to declare
Scientific Advisory Board	No relevant conflicts of interest to declare

Presentation includes discussion of the following off-label use of a drug or medical device:

<N/A>

Introduction



Metrology for health call 2011

- support reliable and efficient exploitation of diagnostic and therapeutic techniques and development of new technologies to improve healthcare
- metrology is the science of measurement





Dutch
Metrology
Institute



Metrological characterization of microvesicles from body fluids as non-invasive diagnostic biomarkers



METAS

Federal Office
of Metrology



Letters of Support



Children's Hospital Boston

universiteit
hasselt

SickKids®

DANA-FARBER
CANCER INSTITUTE



University of Twente
The Netherlands

GENFIT
TOWARDS BETTER MEDICINE

Stago
Education Available
24 Hours a Day 7 Days a Week

Oxford Radcliffe Hospitals NHS
NHS Trust

Maastricht University




University Medical Center
Utrecht

VU medisch centrum




ACTA
academisch centrum tandheelkunde amsterdam



Aim


- develop reliable, comparable and quantitative analysis of microvesicles in biological fluids
 - development of isolation procedures
 - dimensional characterization
 - characterization of the chemical composition, morphology and concentration
 - selection, characterization and distribution of reference materials

Work packages




Work package 1

- development and application of procedures for microvesicle isolation



Work package 2




- dimensional characterization of microvesicles and reference materials
 - free in suspension (nanoparticle tracking analysis, resistive pulse sensing, small angle X-ray scattering)
 - adhered to a surface
 - dried conditions (atomic force microscopy, (transmission) scanning electron microscopy)
 - wet conditions (atomic force microscopy)





BESSYII

Small angle X-ray scattering



Work package 3

- chemical analysis, morphology, and concentration of microvesicles
 - chemical analysis (anomalous small angle x-ray scattering, X-ray fluorescence)
 - cellular origin and type (atomic force microscopy with functionalized tips)
 - morphology
 - dried conditions (atomic force microscopy, transmission electron microscopy)
 - wet conditions (atomic force microscopy)
 - concentration (nanoparticle tracking analysis, resistive pulse sensing)

Work package 4

- development and distribution of traceable reference materials

Reference material	Size (nm)	Concentration (ml ⁻¹)	Density (g/cm ³)	Refractive index @530 nm
Synthetic particles				
• polystyrene beads	30 – 1,000	$1 \times 10^{10} - 1 \times 10^{14}$	1.05	1.599
• silica beads	30 – 1,000	$1 \times 10^{10} - 1 \times 10^{14}$	2.00	1.461
biological particles				
• Intralipid	25 – 700	$\sim 1 \times 10^{14}$	0.93	1.465
• purified vesicles	30 – 1,000	variable	1.13-1.19	not known

- inter metrological laboratory comparison
- inter clinical laboratory comparison

Inter clinical laboratory comparison

- goal
 - validate developed protocols and detection methods in clinical laboratories using traceable reference materials
- distribution
 - September and December 2014
- data
 - collection and analysis in January and February 2015
- results
 - report and peer-reviewed article

Participation

To participate in this SSC survey, please send an e-mail to

r.nieuwland@amc.uva.nl

before October 1st 2012. Please include your

- name and affiliations
- available detection method(s)